

SURGICAL TECHNIQUE









# LLIF Procedure

Lateral Lumbar Interbody Fusion



with MARS™3V & TransContinental®









## Life moves us

At Globus, we move with a sense of urgency to deliver innovations that improve the quality of life for patients with spinal disorders. We are inspired by the needs of these patients and also the needs of the surgeons and health care providers who treat them.

This passion combined with Globus' world class engineering transforms clinical insights into tangible spine care solutions. We are driven to provide the highest quality products to improve the techniques and outcomes of spine surgery so patients can resume their lives as quickly as possible. We extend our reach beyond our world class implants, instrumentation, and service by partnering with researchers and educators to advance the science and knowledge of spine care.


The energy and enthusiasm each of us bring everyday to Globus is palpable. We are constantly in the pursuit of better patient care and understand that speed is critical because life cannot wait.







### Lateral Lumbar Interbody Fusion



During the past 35 years, surgeons have continued to develop the lateral technique. In 1976 Fraser et al. described a muscle splitting approach to the lumbosacral spine. The exposure gained with this retroperitoneal technique offered successful access to the mid-lumbar and lumbosacral spine.3 O'Brien analyzed his many patients that underwent a leftsided, retroperitoneal approach. The benefit of this was avoidance of canal trauma and elimination of motion by using a large graft.<sup>4</sup> McAfee et al. presented the first clinical series investigating the use of the lateral retroperitoneal approach for fusion from L1-L5, reporting an overall lower morbidity comparable to traditional approaches.<sup>5</sup> Building on this, Bertagnoli developed the anterolateral transpsoatic approach for accessing the lumbar disc while avoiding disruption to the posterior elements.6 Today, Globus has refined the lateral approach with the introduction of MARS™3V and TransContinental®.

- 1 Ménard V: Causes de la paraplégie dans le mal de Pott. Rev Orthop: 47-64, 1894.
- 2 Capener N: The evolution of lateral rhacotomy. J Bone Joint Surg Br 36: 173-179, 1954.
- 3 Fraser RD, Gogan WJ. A Modified Muscle-Splitting Approach to the Lumbosacral Spine. Spine 17(8): 943-8, 1992.
- 4 Kozak I. O'Brien I. Simultaneous Combined Anterior and Posterior Fusion. An Independent Analysis of a Treatment for the Disabled Low-Back Pain Patient Spine 15(4): 322-328, 1990.
- 5 McAfee P et al. Minimally Invasive Anterior Retroperitoneal Approach to the Lumbar Spine: Emphasis on the Lateral BAK. Spine 23(13): 1476-1484, 1998.
- 6 Bertagnoli R, Vazquez RJ. The Anterolateral TransPsoactic Approach(ALPA): a new technique for implanting prosthetic disc-nucleus devices. J Spinal Disord Tech 16(4): 398-404, 2003.



### LATERAL LUMBAR INTERBODY FUSION

The MIS lateral approach has been refined with the combination of the MARS™3V retractor and the TransContinental® spacer system.

MARS™3V is a versatile and variable vision system that gives unprecedented control of tissue retraction during the entire procedure.

TransContinental® is a comprehensive spacer system with extensive instrumentation.

### MARS™3V


The MARS™3V Retractor provides the versatility, stability, and precision required to access hard-to-reach disc spaces.

- Independent blade retraction
- Individual blade angulation up to 20°
- Lightweight, aluminum components for radiolucency
- Illumination system fits within blade walls for less bulk

### TransContinental®

The TransContinental® System is ideal for a minimally invasive approach to help preserve patient anatomy.

- Self-distracting leading edge for ease of insertion
- Radiographic positioning markers for implant placement and orientation
- Comprehensive disc preparation instrumentation
- Trial length indicators provide accurate length assessment on fluoroscopy






## **CONTENTS**

Distinguishing Characteristics	2
MARS™3V Retractor Overview	4
Implant Overview	5
Instrument Overview	6
Surgical Technique	
1. Patient Preparation	16
2. Initial Dilator Insertion	18
3. Retractor Insertion	19
4. Disc Preparation	25
5. Interbody Insertion	26
Using the Implant Insertion Tool	27
Final Position	28
TransContinental® Implant Set	30
MARS™3V Retractor Set	32
MARS™ Instrument II Set	34
Anterior Disc Prep I Instrument Set	36
Anterior Disc Prep II Instrument Set	38
Lateral Disc Prep Instrument Set	40
TransContinental® Insertion Instruments Set	42
Lateral Angled Disc Prep Instrument Set	44
TransContinental® Angled Trials Set	46
TransContinental® 16mm Wide Implant Set	48
Bayonetted Disc Prep Instrument Set	49
TransContinental® 22mm Implant Set	50
TransContinental® 10° Lordotic Implant Set	52
TransContinental® Coronal Tapered Implant Set	54
Important Information	56

The Surgical Technique shown is for illustrative purposes only. The technique(s) actually employed in each case always depends on the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient. Additionally, as instruments may occasionally be updated, the instruments depicted in this Surgical Technique may not be exactly the same as the instruments currently available. Please consult with your sales representative or contact Globus directly for more information.

## MARS™3V RETRACTOR OVERVIEW




### **IMPLANT OVERVIEW**


### TransContinental® Spacer


- Self-distracting leading edge
- Large single chamber for bone graft material
- Five axial footprints
- Three sagittal profiles (0°, 6° and 10° Lordotic)
- Heights of 7–17 mm in 2mm increments
- Pyramidal teeth resist migration


### TransContinental® Instrumentation

- MIS compatible system
- Each trial provides precise length measurements
- Holder with sleeve eliminates challenges with implant disengagement










## **INSTRUMENT OVERVIEW**



## Retractor Assembly

Retractor Blades				
Length	Posterior Blade	CC Blade		
40mm	698.450	698.510		
50mm	698.452	698.512		
60mm	698.454	698.514		
70mm	698.456	698.516		
80mm	698.458	698.518		
90mm	698.460	698.520		
100mm	698.462	698.522		
110mm	698.464	698.524		
120mm	698.466	698.526		
130mm	698.468	698.528		
140mm	698.470	698.530		
150mm	698.472	698.532		
160mm	698.474	698.534		
170mm	698.476	698.536		



**Posterior Blade** 




CC Blade




Retractor 3 Blade Frame 698.100



Table Clamp 632.500



### Retractor Instruments



## Light Components





Adapter, ACMI 632.305



Adapter, Olympus 632.307




Adapter, Wolf 632.306



Adapter, Storz 632.308

## Disc Prep Instruments



## Disc Prep Instruments (cont'd)



Thin Rasp, 12x20mm 675.510




Cobb, 7° Up-Angled, 10mm 675.515



Cobb, 7° Up-Angled, 10mm 675.516

Kerrison, 4mm 625.202 Kerrison, 6mm 625.203



Disc Rongeur, 4mm, Up Biting 625.306

Disc Rongeur, 6mm 625.307



## Scrapers



 Height	Part Number
7mm	675.607
9mm	675.609
11mm	675.611
13mm	675.613
15mm	675.615
17mm	675.617

### Curettes



Cobb Elevator, Angled, 18mm 625.102



Ring Curette, 10mm 625.401



Ring Curette, 15mm 625.402

### Curettes (cont'd)



Bone Curette, 7.5 x 11.5mm, Straight 625.407



Bone Curette, 7.5 x 11.5mm, Up-Angled 625.408



Ring Curette, 10mm, Straight 675.518



Ring Curette, 10mm, 7° Up-Angled 675.519



Cup Curette, 6.5 x 9.5mm, Straight 675.525




Cup Curette, 6.5 x 9.5mm, 15° Up-Angled 675.526



Cup Curette,  $6.5 \times 9.5 \text{mm}$ ,  $90^{\circ}$  Down-Angled 675.527

## Handles



### **Insertion Instruments**





Trial, Lordotic

Trials					
Height	Parallel	Lordotic 6°			
5mm	675.006	675.106			
7mm	675.007	675.107			
9mm	675.009	675.109			
11mm	675.011	675.111			
13mm	675.013	675.113			
15mm	675.015	675.115			
17mm	675.017	675.117			

## Insertion Instruments (cont'd)



Slide Hammer, Long 675.004



Implant Insertion Tool 664.500



Insertion Sleeve 675.501



Implant Insertion Tool 664.500, Insertion Sleeve 675.501 (Assembled)



Implant Positioner 675.502

## LLIF SURGICAL TECHNIQUE

## Step 1


### **Patient Preparation**

#### **Patient Positioning**

The patient is placed on a flexible surgical table in a true 90° right lateral decubitus position so that the iliac crest is just over the table break, as shown below.

The patient is then secured to the table at the following locations: (1) Just beneath the iliac crest; (2) Over the thoracic region, just under the shoulder; (3) From the back of the table, over the ankle, and (4) past the knee to the front of the table.

The table should be flexed to open the interval between the 12th rib and iliac crest, and provide direct access to the disc space as shown below.



Patient positioning




Table flexed

#### X-Ray Confirmation

Fluoroscopy is used to ensure that the spine is oriented in a true lateral position. The table should be adjusted so that the C-arm provides true AP images when at 0° and true lateral images at 90°.






Lateral image


AP image

#### **Incision Location**

The operative area is carefully cleaned and the **Incision Locator** is used under fluoroscopy to identify the middle of the disc space to be fused. An access incision mark is then traced on the patient's skin to indicate the position and insertion site for the retractor.



**Using Incision Locator** 




Marking the incision locations

### **Initial Dilator Insertion**

The initial dilator is inserted and AP fluoroscopy is used to confirm proper initial dilator alignment. A K-wire is inserted through the initial dilator into the disc space, in preparation for sequential dilation.

The **Initial Dilator Holder** may be used to steady the dilator and maintain distance from x-ray field while imaging.





Initial dialor position


**Sequential Dilation** 

With the K-wire in place, a series of cannulas are passed over the initial dilator, spreading the tissue to prepare for retractor insertion.

Cannula D should only be used if additional dilation is necessary following Cannula B, C. However, if Cannula D is used, it must be removed prior to insertion of the retractor.

### **Blade Length**

Depth markings on each cannula are used to estimate retractor blade length. The targeted depth indicator will be the first visible marking above the incision site. To ensure sufficient length, 10mm should be added to the blade length reading.



## Step

### **Retractor Insertion**

#### **Retractor Assembly**

Select the appropriate blades and insert into each of the three blade mounts of the Retractor 3 Blade Frame.

The **Posterior Blade** is inserted into the posterior blade mount.

The CC Blades (cephalad-caudal) are inserted into the cephalad and caudal blade mounts.







CC Blade



Ensure that the blades are properly seated into the retractor at each of the three positions. Secure the blades using the **Hook and Latch Driver**. Position the driver on the latch and rotate 90° clockwise to lock the blade in place.

The blades can be changed intraoperatively when a different blade length is required. The blades have angled holes that accept the driver. Insert the hook and tighten down the white sleeve to hold the blade securely. This provides a secure connection to remove the blade.



Using the Hook and Latch Driver to secure blades



Using the Hook and Latch Driver to change blades

### Retractor Insertion (cont'd)

#### **Retractor Positioning**

Ensure that the retractor is in the fully closed position and the blades are securely attached to the frame. The access incision should allow the blades to retract and angulate.

Slide the retractor over Cannula C and apply gentle downward pressure on the frame.

Before removing the cannulas, angulate all three blades to one full turn of the silver knobs. Retract all three blades to two clicks using the gold knobs, starting with the blade closest to the iliac crest. Angulating and retracting the blades in this manner will help prevent tissue creep as the cannulas are removed.

Once the retractor has been securely positioned and the Articulating Arm Assembly tightened, remove the cannulas and verify position of the retractor before removing the k-wire.

Use AP fluoroscopy to verify the correct positioning on the spine and to confirm that retractor blades are parallel with the disc space.


#### **Table Arm Attachment**

Attach the **Table Clamp** over the drape and onto the bed rail attachment. Insert the Articulating Arm Assembly into the clamp and secure.


The opposite end of the assembly arm is then attached to the retractor. There are two options for attachment positions on the retractor, as shown at right.

Attaching the arm assembly to point A maintains retractor position relative to the posterior blade position and translates the cephalad and caudad blades anteriorly when the retractor is opened.

Attaching the arm assembly to point B maintains the retractor position relative to the center of the frame and retracts all three blades when the retractor is opened.



Positioning retractor



Articulating Arm attachment locations

#### Table Arm Attachment (cont'd)

Insert the Articulating Arm Assembly into the desired attachment position and tighten the thumb screw using the 10mm **Socket Driver**. Position the arm and lock in place by tightening the black star handle on the arm assembly.

Minimal torque is required to tighten the thumb screw with the driver.

Manipulation of the retractor can be achieved with the Frame Handle that fits over arm attachment point B.



Using a Frame Handle

### **Light Cable Insertion**

The MIS Illumination System has two light cables which can be inserted into the CC Blades.

The light cables should be inserted through the blade to a depth providing optimal visibility.

The **Fiber Optic Cord** attaches to the light source used for head lamps or endoscopes. The adapters accommodate an ACMI, Olympus, Storz, or Wolf light source.

Note: The Fiber Optic Cord is reusable and must not be discarded.



### Retractor Insertion (cont'd)

### Cephalad-Caudad Blade Anchoring

For additional retractor stability, **Docking Pins** can be inserted into the vertebral body through the CC Blades to increase retractor stability when expanding blades for greater exposure as shown below.

Select the appropriate size pin and insert into the **Docking Pin Sleeve**. The pin assembly can slide down the T-slot on either side of the blade. The Docking Pin **Tool** has a hex feature that mates with the head of the pin. Rotate the tool clockwise to engage pin threads into the bone.

To remove the pins, re-engage the hex of the tool into the pin head and rotate counterclockwise to disengage from the bone.



Docking pin assembly



Docking pin assembly loaded



#### **Blade Expansion**

Each blade may be independently expanded using the 10mm Socket Driver to rotate the respective gold hex nut in the direction indicated by the arrow.

Each blade may be angled up to  $20^\circ$  using the driver. Place the driver onto the silver hex nut and rotate the instrument clockwise, allowing the blade to tilt to the desired position.




### Retractor Insertion (cont'd)

#### Blade Expansion (cont'd)

#### Widening Shim

While expanding the CC Blades, Widening Shims can be used to help prevent soft tissue creep between blades.

The **Shim Tool, CC**, is used to slide the Widening Shim down the T-slot on either side of the blades, for an additional 22mm of blade width.



#### **Lengthening Shim**


While angulating the CC Blades, Lengthening Shims can be used to help prevent soft tissue creep between blades by increasing the length of the blade to maintain bone contact.

The Shim Tool, CC, is used to push the Lengthening Shim down the right T-slot of the blades.



#### Shim Removal

Both Widening and Lengthening Shims have angled holes to accept the hook on the Shim Tool, CC. Insert the hook and pull upward to remove the shim.






## **Disc Preparation**

### **Annulotomy**

The Bayonneted Annulotomy Knife is used to create a window centered in the anterior half of the annulus large enough for graft insertion.

#### **Contralateral Annulus Release**


A **Cobb Elevator** is passed along both endplates through the disc space, far enough to provide release of the contralateral annulus. This allows for height restoration upon insertion of the implant.



Using the Cobb Elevator

### **Disc Space Preparation**

Leaving the posterior annulus intact, remove the intervertebral disc and osteophytes as needed. The Disc Box Cutter, Rotary Cutter, Disc Rongeurs, Kerrisons, Curettes, Scrapers and Rasps are provided for disc removal and endplate preparation, as shown at right.



Using the Disc Rongeur


## Step

### **Interbody Insertion**

### **Implant Sizing**

To determine the appropriate spacer size for the desired segment, first insert the smallest **Trial** into the disc space, moving to larger trials as needed. Determine which trial best fits the prepared disc space. A secure fit is desirable in order to maintain disc height and stabilize the segment. Use AP fluoroscopy to confirm that the implant is centered and lateral fluoroscopy to ensure that the implant is in the appropriate AP position.

The trial's anterior side is labeled "Anterior" and matches the profile of the TransContinental® Spacer.






**Implant Sizing** (55mm disc space is shown)


### **Implant Insertion**

Select the appropriately sized TransContinental® Spacer and gently insert into the intervertebral space using the Implant Holder and Tamp. AP fluoroscopy should be used to facilitate implant placement.

Once the position is confirmed, the implant is released from the holder.


The TransContinental® Spacer is intended to be used with supplemental fixation.





#### Closure

Once the procedure is completed, use the 10mm Socket Driver to angle all blades to the 0° starting position. Return all blades to the closed position by compressing the three release tabs on the back and sides of the retractor.



## **Final Position**




AP view (radiographic markers shown)



Sagittal view (radiographic markers shown)





Axial view



Radiographic markers positions

## TransContinental® IMPLANT SET








## TransContinental® Implant Set 975.916

Part No.	Description	Qty	Part No.	Description	Qty
375.007	TransContinental® Spacer, Small, 0°, 7mm	2	375.513	TransContinental® Spacer, Large, 6°, 13mm	2
375.009	TransContinental® Spacer, Small, 0°, 9mm	2	375.515	TransContinental® Spacer, Large, 6°, 15mm	2
375.011	TransContinental® Spacer, Small, 0°, 11mm	2	375.517	TransContinental® Spacer, Large, 6°, 17mm	1
375.013	TransContinental® Spacer, Small, 0°, 13mm	2	375.607	TransContinental® Spacer, X-Large, 0°, 7mm	2
375.015	TransContinental® Spacer, Small, 0°, 15mm	2	375.609	TransContinental® Spacer, X-Large, 0°, 9mm	2
375.017	TransContinental® Spacer, Small, 0°, 17mm	1	375.611	TransContinental® Spacer, X-Large, 0°, 11mm	1 2
375.107	TransContinental® Spacer, Small, 6°, 7mm	2	375.613	TransContinental® Spacer, X-Large, 0°, 13mm	1 2
375.109	TransContinental® Spacer, Small, 6°, 9mm	2	375.615	TransContinental® Spacer, X-Large, 0°, 15mm	1 2
375.111	TransContinental® Spacer, Small, 6°, 11mm	2	375.617	TransContinental® Spacer, X-Large, 0°, 17mm	n 1
375.113	TransContinental® Spacer, Small, 6°, 13mm	2	375.707	TransContinental® Spacer, X-Large, 6°, 7mm	2
375.115	TransContinental® Spacer, Small, 6°, 15mm	2	375.709	TransContinental® Spacer, X-Large, 6°, 9mm	2
375.117	TransContinental® Spacer, Small, 6°, 17mm	1	375.711	TransContinental® Spacer, X-Large, 6°, 11mm	1 2
375.207	TransContinental® Spacer, Medium, 0°, 7mm	2	375.713	TransContinental® Spacer, X-Large, 6°, 13mm	1 2
375.209	TransContinental® Spacer, Medium, 0°, 9mm	2	375.715	TransContinental® Spacer, X-Large, 6°, 15mm	1 2
375.211	TransContinental® Spacer, Medium, 0°, 11mm	2	375.717	TransContinental® Spacer, X-Large, 6°, 17mm	n 1
375.213	TransContinental® Spacer, Medium, 0°, 13mm	2	375.807	TransContinental® Spacer, X-Small, 0°, 7mm	2
375.215	TransContinental® Spacer, Medium, 0°, 15mm	2	375.809	TransContinental® Spacer, X-Small, 0°, 9mm	2
375.217	TransContinental® Spacer, Medium, 0°, 17mm	1	375.811	TransContinental® Spacer, X-Small, 0°, 11mm	2
375.307	TransContinental® Spacer, Medium, 6°, 7mm	2	375.813	TransContinental® Spacer, X-Small, 0°, 13mm	2
375.309	TransContinental® Spacer, Medium, 6°, 9mm	2	375.815	TransContinental® Spacer, X-Small, 0°, 15mm	2
375.311	TransContinental® Spacer, Medium, 6°, 11mm	2	375.817	TransContinental® Spacer, X-Small, 0°, 17mm	1
375.313	TransContinental® Spacer, Medium, 6°, 13mm	2	375.907	TransContinental® Spacer, X-Small, 6°, 7mm	2
375.315	TransContinental® Spacer, Medium, 6°, 15mm	2	375.909	TransContinental® Spacer, X-Small, 6°, 9mm	2
375.317	TransContinental® Spacer, Medium, 6°, 17mm	1	375.911	TransContinental® Spacer, X-Small, 6°, 11mm	. 2
375.407	TransContinental® Spacer, Large, 0°, 7mm	2	375.913	TransContinental® Spacer, X-Small, 6°, 13mm	. 2
375.409	TransContinental® Spacer, Large, 0°, 9mm	2	375.915	TransContinental® Spacer, X-Small, 6°, 15mm	. 2
375.411	TransContinental® Spacer, Large, 0°, 11mm	2	375.917	TransContinental® Spacer, X-Small, 6°, 17mm	ı <b>1</b>
375.413	TransContinental® Spacer, Large, 0°, 13mm	2	975.011	TransContinental® Graphic Case – Implants	
375.415	TransContinental® Spacer, Large, 0°, 15mm	2			
375.417	TransContinental® Spacer, Large, 0°, 17mm	1			
375.507	TransContinental® Spacer, Large, 6°, 7mm	2			
375.509	TransContinental® Spacer, Large, 6°, 9mm	2			
375.511	TransContinental® Spacer, Large, 6°, 11mm	2			

## MARS™3V RETRACTOR SET









### MARS<sup>™</sup>3V Retractor Set 998.901

	Instrume	nts	Qty	Retracto	r Blades	Qty
1	623.003	K-Wire Gripper	1	698.476	Blade, Posterior, 170mm	2
2	698.100	Retractor 3 Blade Frame	1	698.510	Blade, CC, 40mm	2
3	632.102	Retractor 2 Blade Frame	1	698.512	Blade, CC, 50mm	2
4	632.150	10mm Socket Driver	1	698.514	Blade, CC, 60mm	2
5	698.250	Hook and Latch Driver	1	698.516	Blade, CC, 70mm	2
6	675.403	Flouro Modulator	1	698.518	Blade, CC, 80mm	2
7	675.404	Incision Locator	1	698.520	Blade, CC, 90mm	2
8	675.513	8" Suction	1	698.522	Blade, CC, 100mm	2
9	675.800	Radiolucent Initial Dilator Holder	1	698.524	Blade, CC, 110mm	2
10	698.205	Cannula A	1	698.526	Blade, CC, 120mm	2
11	698.210	Cannula B	1	698.528	Blade, CC, 130mm	2
12	698.215	Cannula C	1	698.530	Blade, CC, 140mm	2
13	698.220	Cannula D	1	698.532	Blade, CC, 150mm	2
14	698.230	Frame Handle	1	698.534	Blade, CC, 160mm	2
15	698.240	Shim Tool, CC	1	698.536	Blade, CC, 170mm	2
16	698.260	Docking Pin Tool	1	Disposal	and land	
17	698.330	Disc Shim Tool	1	•	Bipolar Forceps, 10" Bayo, 1.0mm Tip	1
18	698.350	Docking Pin Sleeve	4		MARS™3V Disposable Kit	1
(19)	Retractor	r Rlados			Lengthening Shim	2
	698.450	Blade, Posterior, 40mm	2		Widening Shim	2
	698.452	Blade, Posterior, 50mm	2		Docking Pin, 10mm	2
	698.454	Blade, Posterior, 60mm	2		Docking Pin, 20mm	2
	698.456	Blade, Posterior, 70mm	2		Disc Shim, Aluminum	1
	698.458	Blade, Posterior, 80mm	2		Disc Shim, Stainless Steel	•
	698.460	Blade, Posterior, 90mm	2	070.3203	Disc Stiff, Staffiess Steel	
	698.462	Blade, Posterior, 100mm	2			
	698.464	Blade, Posterior, 110mm	2			
	698.466	Blade, Posterior, 120mm	2			
	698.468	Blade, Posterior, 130mm	2			
	698.470	Blade, Posterior, 140mm	2			
	698.472	Blade, Posterior, 150mm	2			
	698.474	Blade, Posterior, 160mm	2			
	U)U.T/T	Diage, Fosterior, Footiliti	_			

## MARS<sup>™</sup> INSTRUMENT II SET








# MARS<sup>™</sup> Instrument II Set 932.902

	Instruments			
1	632.300	Fiber-Optic Cord	1	
2	632.305	Adapter, ACMI	1	
3	632.306	Adapter, Wolf	1	
4	632.307	Adapter, Olympus	1	
5	632.308	Adapter, Storz	1	
6	632.390	Port Lock, 19mm	1	
7	632.391	Port Lock, 22mm	1	
8	632.392	Port Lock, 26mm	1	
9	632.401	2mm Cannula	1	
10	632.402	5mm Cannula	1	
11	632.403	8mm Cannula	1	
12	632.404	12mm Cannula	1	
13	632.405	15mm Cannula	1	
14	632.406	18mm Cannula	1	
15	632.407	22mm Cannula	1	
16	632.408	26mm Port Mount	1	
17	632.409	26mm Port Positioner	1	
18	632.410	22mm Port Mount	1	
19	632.411	22mm Port Positioner	1	
20	632.412	19mm Port Mount	1	
21	632.413	19mm Port Positioner	1	
22	632.500	Table Clamp	1	
23	632.750	Articulating Arm Assembly	1	
	932.002	MARS™ Instrument II Graphic Case		

# ANTERIOR DISC PREP I INSTRUMENT SET










# Anterior Disc Prep I Instrument Set 925.901

	Instruments				
1	625.201	Kerrison, 2mm	1		
2	625.202	Kerrison, 4mm	1		
3	625.203	Kerrison, 6mm	1		
4	625.301	Bone Rongeur, Doube Acting, 8mm	1		
5	625.302	Bone Rongeur, Double Acting, 12mm	1		
6	625.303	Disc Rongeur, 2mm	1		
7	625.304	Disc Rongeur, 2mm, Up Biting	1		
8	625.305	Disc Rongeur, 4mm	1		
9	625.306	Disc Rongeur, 4mm, Up Biting	1		
10	625.307	Disc Rongeur, 6mm	1		
11	625.609	ALIF Trial Sizer, 9mm	1		
12	625.611	ALIF Trial Sizer, 11mm	1		
13	625.613	ALIF Trial Sizer, 13mm	1		
14	625.615	ALIF Trial Sizer, 15mm	1		
15	625.617	ALIF Trial Sizer, 17mm	1		
16	625.619	ALIF Trial Sizer, 19mm	1		
17	625.621	ALIF Trial Sizer, 21mm	1		
18	625.709	Rotary Scraper, 9mm	1		
19	625.711	Rotary Scraper, 11mm	1		
20	625.713	Rotary Scraper, 13mm	1		
21	625.715	Rotary Scraper, 15mm	1		
22	625.717	Rotary Scraper, 17mm	1		
23	625.801	Vein Retractor	1		
24	625.804	T-Handle with Impaction Cap, Long	2		
25	625.805	Parallel Distractor	1		
26	625.806	Knife Handle	1		
27	625.811	Long Penfield	1		
28	675.005	T-Handle with Impaction Cap	2		
	925.101	Graphic Case			

# ANTERIOR DISC PREP II INSTRUMENT SET










# Anterior Disc Prep II Instrument Set 925.902

	Instruments					
1	625.101	Cobb Elevator, 18mm	1			
2	625.102	Cobb Elevator, Angled, 18mm	1			
3	625.103	Cobb Elevator, 23mm	1			
4	625.104	Cobb Elevator, Angled, 23mm	1			
5	625.401	Ring Curette, 10mm	1			
6	625.402	Ring Curette, 15mm	1			
7	625.403	Bone Curette, 3.5mm x 5.5mm, Straight	1			
8	625.404	Bone Curette, 3.5mm x 5.5mm, Up-Angled	1			
9	625.405	Bone Curette, 5.5mm x 8.5mm, Straight	1			
10	625.406	Bone Curette, 5.5mm x 8.5mm, Up-Angled	1			
11	625.407	Bone Curette, 7.5mm x 11.5mm, Straight	1			
12	625.408	Bone Curette, 7.5mm x 11.5mm, Up-Angle	d 1			
13	625.501	Dual Rasp	1			
14	625.502	Angled Rasp	1			
15	625.803	Osteotome, 16mm x 20mm	1			
	925.102	Graphic Case II				
	Addition	ally Available				
	625.409	Bone Curette, 9.5mm x 14.5mm, Straight				
	625.410	Bone Curette, 9.5mm x 14.5mm, Up-Angle	d			
	625.411	Bone Curette, 11.5mm x 17.5mm, Straight				
	625.412	Bone Curette, 11.5mm x 17.5mm, Up-Ang	led			
	625.413	Bone Curette, 13.5mm x 20.5mm, Straight				
	625.414	Bone Curette, 13.5mm x 20.5mm, Up-Ang	led			

# LATERAL DISC PREP INSTRUMENT SET







# Lateral Disc Prep Instrument Set 975.914


	Instrume	ents	Qty		
1	675.002	Slap Hammer Adaptor	1		
2	675.005 T-Handle with Impaction Cap				
3	675.201	Quick Connect Guide	2		
4	675.503	Large Cobb Elevator	1		
5	675.504	Small Cobb Elevator	1		
6	675.505	Rasp	1		
7	675.506	Disc Rotary Cutter	1		
8	675.507	Box Cutter	1		
9	675.510	Thin Rasp, 12x20mm	1		
10	675.515	Cobb, 10mm, 7°, Up-Angle	1		
11	675.516	Cobb, 20mm, 7°, Up-Angle	1		
12	675.518	Ring Curette, 10mm, Straight	1		
13	675.519	Ring Curette, 10mm, 7°, Up-angle Tip	1		
14	675.525	Cup Curette, 6.5x9.5mm, Straight	1		
15	675.526	Cup Curette, 6.5x9.5mm, 15°, Up-angle	1		
16	675.527	Cup Curette, 6.5x9.5mm, 90°, Down-angle	1		
17	675.607	Scaper, 7mm	1		
18	675.609	Scaper, 9mm	1		
19	675.611	Scaper, 11mm	1		
20	675.613	Scaper, 13mm	1		
21	675.615	Scaper, 15mm	1		
22	675.617	Scaper, 17mm	1		
	975.008	TransContinental® Disc Preparation Graphic Case			


# Additionally Available


675.170S Bipolar Forceps Bayonetted, Straight

675.171S Bipolar Forceps Bayonetted, Angled

# TransContinental® INSERTION INSTRUMENTS SET










# TransContinental® Insertion Instruments Set 975.915

	Instruments					
1	664.500	PATRIOT® CONTINENTAL® Holder/Inserter	2			
2	675.004	Long Throw Slide Hammer	1			
3	675.005	T-Handle with Impaction Cap	2			
4	675.006	TransContinental® Trial, Parallel, 5mm	1			
5	675.007	TransContinental® Trial, Parallel, 7mm	1			
6	675.009	TransContinental® Trial, Parallel, 9mm	1			
7	675.011	TransContinental® Trial, Parallel, 11mm	1			
8	675.013	TransContinental® Trial, Parallel, 13mm	1			
9	675.015	TransContinental® Trial, Parallel, 15mm	1			
10	675.017	TransContinental® Trial, Parallel, 17mm	1			
11	675.106	TransContinental® Trial, Lordotic, 5mm	1			
12	675.107	TransContinental® Trial, Lordotic, 7mm	1			
13	675.109	TransContinental® Trial, Lordotic, 9mm	1			
14	675.111	TransContinental® Trial, Lordotic, 11mm	1			
15	675.113	TransContinental® Trial, Lordotic, 13mm	1			
16	675.115	TransContinental® Trial, Lordotic, 15mm	1			
17	675.117	TransContinental® Trial, Lordotic, 17mm	1			
18	675.501	Insertion Sleeve	2			
19	675.502	Tamp	1			
20	679.010	L-Handle	1			
	975.007	TransContinental® Graphic Case – Insertion				

# LATERAL ANGLED DISC PREP INSTRUMENT SET









# Lateral Angled Disc Prep Instrument Set 975.940

	Instruments				
1	675.005	T-Handle with Impaction Cap	2		
2	675.330	Cup, Curette, Angled, 6.5x9.5mm, Down	1		
3	675.331	Cup, Curette, Angled, 6.5x9.5mm, Up	1		
4	675.334	Cobb, Angled, 10mm, Down	1		
5	675.335	Cobb, Angled, 10mm, Up	1		
6	675.336	Cobb, Angled, 20mm, Down	1		
7	675.337	Cobb, Angled, 20mm, Up	1		
8	675.338	Double Sided Rasp, Angled	1		
9	675.339	Double Sided Scraper, Angled	1		
10	675.340	Rongeur, Angled, Left, 4mm	1		
11	675.341	Rongeur, Angled Right, 4mm	1		
12	675.342	Rongeur, Angled Left, 6mm	1		
13	675.343	Rongeur, Angled Right, 6mm	1		
14	675.344	Box Cutter, Angled, 6x18mm	1		
15	675.345	Box Cutter, Angled, 8x18mm	1		
16	675.346	Box Cutter, Angled, 10x18mm	1		
17	675.347	Tamp, Angled	1		
18	675.980	Angled Inserter	2		
	975.040	Lateral Angled Disc Prep Graphic Case			

# TransContinental® Angled Trials Set





# TransContinental® Angled Trials Set 975.941

Part No.	Description	Qty
675.273	18mm Trial, Angled, Parallel, 5mm	1
675.275	18mm Trial, Angled, Parallel, 7mm	1
675.277	18mm Trial, Angled, Parallel, 9mm	1
675.279	18mm Trial, Angled, Parallel, 11mm	1
675.281	18mm Trial, Angled, Parallel, 13mm	1
675.473	22mm Trial, Angled, Parallel, 5mm	1
675.475	22mm Trial, Angled, Parallel, 7mm	1
675.477	22mm Trial, Angled, Parallel, 9mm	1
675.479	22mm Trial, Angled, Parallel, 11mm	1
675.481	22mm Trial, Angled, Parallel, 13mm	1
675.543	18mm Trial, Angled, Lordotic, 5mm	1
675.545	18mm Trial, Angled, Lordotic, 7mm	1
675.547	18mm Trial, Angled, Lordotic, 9mm	1
675.549	18mm Trial, Angled, Lordotic, 11mm	1
675.551	18mm Trial, Angled, Lordotic, 13mm	1
675.573	22mm Trial, Angled, 6° Lordotic, 5mm	1
675.575	22mm Trial, Angled, 6° Lordotic, 7mm	1
675.577	22mm Trial, Angled, 6° Lordotic, 9mm	1
675.579	22mm Trial, Angled, 6° Lordotic, 11mm	1
675.581	22mm Trial, Angled, 6° Lordotic, 13mm	1
975.041	TransContinental® Angled Trials Graphic Case	
Additiona	ally Available	
675.083	18mm Trial, Angled, Right Lordotic, 5mm	
675.085	18mm Trial, Angled, Right Lordotic, 7mm	
675.087	18mm Trial, Angled, Right Lordotic, 9mm	
675.089	18mm Trial, Angled, Right Lordotic, 11mm	
675.091	18mm Trial, Angled, Right Lordotic, 13mm	
675.673	22mm Trial, Angled, Right Lordotic, 5mm	
675.675	22mm Trial, Angled, Right Lordotic, 7mm	
675.677	22mm Trial, Angled, Right Lordotic, 9mm	
675.679	22mm Trial, Angled, Right Lordotic, 11mm	
675.681	22mm Trial, Angled, Right Lordotic, 13mm	

# TransContinental® 16mm Wide Implant Set 975.920

Part No.	Description	Qty
375.060	TransContinental® Spacer, 0°, 16x25mm, 7mm	2
375.061	TransContinental® Spacer, 0°, 16x25mm, 9mm	2
375.062	TransContinental® Spacer, 0°, 16x25mm, 11mm	2
375.063	TransContinental® Spacer, 0°, 16x25mm, 13mm	2
375.160	TransContinental® Spacer, 0°, 16x30mm, 7mm	2
375.161	TransContinental® Spacer, 0°, 16x30mm, 9mm	2
375.162	TransContinental® Spacer, 0°, 16x30mm, 11mm	2
375.163	TransContinental® Spacer, 0°, 16x30mm, 13mm	2
375.260	TransContinental® Spacer, 0°, 16x35mm, 7mm	2
375.261	TransContinental® Spacer, 0°, 16x35mm, 9mm	2
375.262	TransContinental® Spacer, 0°, 16x35mm, 11mm	2
375.263	TransContinental® Spacer, 0°, 16x35mm, 13mm	2
375.360	TransContinental® Spacer, 0°, 16x40mm, 7mm	2
375.361	TransContinental® Spacer, 0°, 16x40mm, 9mm	2
375.362	TransContinental® Spacer, 0°, 16x40mm, 11mm	2
375.363	TransContinental® Spacer, 0°, 16x40mm, 13mm	2
375.460	TransContinental® Spacer, 0°, 16x20mm, 7mm	2
375.461	TransContinental® Spacer, 0°, 16x20mm, 9mm	2
375.462	TransContinental® Spacer, 0°, 16x20mm, 11mm	2
375.463	TransContinental® Spacer, 0°, 16x20mm, 13mm	2
675.059	16mm Trial 0° Lordotic, 5mm	1
675.060	16mm Trial 0° Lordotic, 7mm	1
675.061	16mm Trial 0° Lordotic, 9mm	1
675.062	16mm Trial 0° Lordotic, 11mm	1
675.063	16mm Trial 0° Lordotic, 13mm	1
675.930	Insertion Sleeve	2
675.940	Holder Inserter	2
975.020	TransContinental® Graphic Case - Implant, 16mm Wide	

# Bayonetted Disc Prep Instrument Set 975.923

Instruments		
675.128	Dual Convex Rasp, Bayonetted	1
675.131	Bone Curette, Bayonetted, Small	1
675.132	Bone Curette, Bayonetted, Up-Angled, Small	1
675.133	Bone Curette, Bayonetted, 90°, Small	1
675.139	Bone Curette, Bayonetted, Large	1
675.140	Bone Curette, Bayonetted, Up-Angled, Large	1
675.141	Bone Curette, Bayonetted, 90°, Large	1
675.148	Ring Curette, Bayonetted, Up Angled	1
675.151	Box Curette, Bayonetted, Up Angled, Wide	1
675.154	Box Curette, Bayonetted, Up Angled, Narrow	1
675.162	Nerve Retractor	1
675.172	Penfield, Bayonetted, Push	1
675.173	Penfield, Bayonetted, Pull	1
675.182	Thin Rasp, Bayonetted, Straight	1
675.183	Thin Rasp, Bayonetted, Angled	1
675.185	Rake, Bayonetted	1
675.186	Nerve Hook, Bayonetted	1
675.305	Kerrison, Bayonetted, 3mm	1
675.306	Kerrison, Bayonetted, 5mm	1
675.316	Disc Rongeur, Bayonetted, 2mm, Straight	1
675.317	Disc Rongeur, Bayonetted, 2mm, Up	1
675.318	Disc Rongeur, Bayo, 4mm, Straight	1
675.319	Disc Rongeur, Bayonetted, 4mm, Up	1
675.320	Disc Rongeur, Bayonetted, 6mm, Straight	1
675.321	Disc Rongeur, Bayonetted, 6mm, Up	1
975.023	Bayonetted Lateral Instruments Graphic Case	

# TransContinental® 22mm Implant Set 975.922

Part No.	Description	Qty	Part No.	Description	Qty
375.045	TransContinental® Spacer, 22mm Wide, Small, 0°, 7mm	1	375.357	TransContinental® Spacer, 22mm Wide, Medium, 6°, 17mm	1
375.049	TransContinental® Spacer, 22mm Wide, Small, 0°, 9mm	2	375.445	TransContinental® Spacer, 22mm Wide, Large, 0°, 7mm	1
375.051	TransContinental® Spacer, 22mm Wide, Small, 0°, 11mm	2	375.449	TransContinental® Spacer, 22mm Wide, Large, 0°, 9mm	2
375.053	TransContinental® Spacer, 22mm Wide, Small, 0°, 13mm	2	375.451	TransContinental® Spacer, 22mm Wide, Large, 0°, 11mm	2
375.055	TransContinental® Spacer, 22mm Wide, Small, 0°, 15mm	2	375.453	TransContinental® Spacer, 22mm Wide, Large, 0°, 13mm	2
375.057	TransContinental® Spacer, 22mm Wide, Small, 0°, 17mm	1	375.455	TransContinental® Spacer, 22mm Wide, Large, 0°, 15mm	2
375.145	TransContinental® Spacer, 22mm Wide, Small, 6°, 7mm	1	375.457	TransContinental® Spacer, 22mm Wide, Large, 0°, 17mm	1
375.149	TransContinental® Spacer, 22mm Wide, Small, 6°, 9mm	2	375.545	TransContinental® Spacer, 22mm Wide, Large, 6°, 7mm	1
375.151	TransContinental® Spacer, 22mm Wide, Small, 6°, 11mm	2	375.549	TransContinental® Spacer, 22mm Wide, Large, 6°, 9mm	2
375.153	TransContinental® Spacer, 22mm Wide, Small, 6°, 13mm	2	375.551	TransContinental® Spacer, 22mm Wide, Large, 6°, 11mm	2
375.155	TransContinental® Spacer, 22mm Wide, Small, 6°, 15mm	2	375.553	TransContinental® Spacer, 22mm Wide, Large, 6°, 13mm	2
375.157	TransContinental® Spacer, 22mm Wide, Small, 6°, 17mm	1	375.555	TransContinental® Spacer, 22mm Wide, Large, 6°, 15mm	2
375.245	TransContinental® Spacer, 22mm Wide, Medium, 0°, 7mm	1	375.557	TransContinental® Spacer, 22mm Wide, Large, 6°, 17mm	1
375.249	TransContinental® Spacer, 22mm Wide, Medium, 0°, 9mm	2	375.645	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 7mm	1
375.251	TransContinental® Spacer, 22mm Wide, Medium, 0°, 11mm	2	375.649	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 9mm	2
375.253	TransContinental® Spacer, 22mm Wide, Medium, 0°, 13mm	2	375.651	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 11mm	2
375.255	TransContinental® Spacer, 22mm Wide, Medium, 0°, 15mm	2	375.653	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 13mm	2
375.257	TransContinental® Spacer, 22mm Wide, Medium, deg, 17mm	1	375.655	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 15mm	2
375.345	TransContinental® Spacer, 22mm Wide, Medium, 6°, 7mm	1	375.657	TransContinental® Spacer, 22mm Wide, X-Large, 0°, 17mm	1
375.349	TransContinental® Spacer, 22mm Wide, Medium, 6°, 9mm	2	375.745	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 7mm	1
375.351	TransContinental® Spacer, 22mm Wide, Medium, 6°, 11mm	2	375.749	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 9mm	2
375.353	TransContinental® Spacer, 22mm Wide, Medium, 6°, 13mm	2	375.751	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 11mm	2
375.355	TransContinental® Spacer, 22mm Wide, Medium, 6°, 15mm	2	375.753	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 13mm	2

# TransContinental® 22mm Implant Set 975.922 (cont'd)

Part No.	Description	Qty	Part No.	Description	Qty
375.755	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 15mm	2	675.453	TransContinental® 22mm Trial, 6°, Lordotic, 13mm	1
375.757	TransContinental® Spacer, 22mm Wide, X-Large, 6°, 17mm	1	675.455	TransContinental® 22mm Trial, 6°, Lordotic, 15mm	1
375.845	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 7mm	1	675.457	TransContinental® 22mm Trial, 6°, Lordotic, 17mm	1
375.849	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 9mm	1	675.522	Insertion Sleeve, 22mm	2
375.851	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 11mm	1	975.022	TransContinental® 22mm Wide Implant Graphic Case	
375.853	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 13mm	1	Additiona	ally Available	
375.855	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 15mm	1	375.857	TransContinental® Spacer, 22mm Wide, X-Small, 0°, 17mm	
375.945	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 7mm	1	375.957	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 17mm	
375.949	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 9mm	1			
375.951	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 11mm	1			
375.953	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 13mm	1			
375.955	TransContinental® Spacer, 22mm Wide, X-Small, 6°, 15mm	1			
675.043	TransContinental® 22mm Trial, Parallel, 5mm	1			
675.045	TransContinental® 22mm Trial, Parallel, 7mm	1			
675.049	TransContinental® 22mm Trial, Parallel, 9mm	1			
675.051	TransContinental® 22mm Trial, Parallel, 11mm	1			
675.053	TransContinental® 22mm Trial, Parallel, 13mm	1			
675.055	TransContinental® 22mm Trial, Parallel, 15mm	1			
675.057	TransContinental® 22mm Trial, Parallel, 17mm	1			
675.443	TransContinental® 22mm Trial, 6° Lordotic, 5mm	1			
675.445	TransContinental® 22mm Trial, 6° Lordotic, 7mm	1			
675.449	TransContinental® 22mm Trial, 6°, Lordotic, 9mm	1			
675.451	TransContinental® 22mm Trial, 6°, Lordotic, 11mm	1			

# TransContinental® 10° Lordotic Implant Set 975.924

Part No.	Description	Qty	Part No.	Description	Qty
375.067	TransContinental® Spacer, 18mm Wide, Small, 10°, 7mm	2	375.377	TransContinental® Spacer, 22mm Wide, Medium, 10°, 17mm	1
375.069	TransContinental® Spacer, 18mm Wide, Small, 10°, 9mm	2	375.467	TransContinental® Spacer, 18mm Wide, Large, 10°, 7mm	2
375.071	TransContinental® Spacer, 18mm Wide, Small, 10°, 11mm	2	375.469	TransContinental® Spacer, 18mm Wide, Large, 10°, 9mm	2
375.073	TransContinental® Spacer, 18mm Wide, Small, 10°, 13mm	2	375.471	TransContinental® Spacer, 18mm Wide, Large, 10°, 11mm	2
375.075	TransContinental® Spacer, 18mm Wide, Small, 10°, 15mm	2	375.473	TransContinental® Spacer, 18mm Wide, Large, 10°, 13mm	2
375.077	TransContinental® Spacer, 18mm Wide, Small, 10°, 17mm	1	375.475	TransContinental Spacer, 18mm Wide, Large, 10°, 15mm	2
375.167	TransContinental® Spacer, 22mm Wide, Small, 10°, 7mm	2	375.477	TransContinental® Spacer, 18mm Wide, Large, 10°, 17mm	1
375.169	TransContinental® Spacer, 22mm Wide, Small, 10°, 9mm	2	375.567	TransContinental® Spacer, 22mm Wide, Large, 10°, 7mm	2
375.171	TransContinental® Spacer, 22mm Wide, Small, 10°, 11mm	2	375.569	TransContinental® Spacer, 22mm Wide, Large, 10°, 9mm	2
375.173	TransContinental® Spacer, 22mm Wide, Small, 10°, 13mm	2	375.571	TransContinental® Spacer, 22mm Wide, Large, 10°, 11mm	2
375.175	TransContinental® Spacer, 22mm Wide, Small, 10°, 15mm	2	375.573	TransContinental® Spacer, 22mm Wide, Large, 10°, 13mm	2
375.177	TransContinental® Spacer, 22mm Wide, Small, 10°, 17mm	1	375.575	TransContinental® Spacer, 22mm Wide, Large, 10°, 15mm	2
375.267	TransContinental® Spacer, 18mm Wide, Medium, 10°, 7mm	2	375.577	TransContinental® Spacer, 22mm Wide, Large, 10°, 17mm	1
375.269	TransContinental® Spacer, 18mm Wide, Medium, 10°, 9mm	2	375.667	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 7mm	2
375.271	TransContinental® Spacer, 18mm Wide, Medium, 10°, 11mm	2	375.669	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 9mm	2
375.273	TransContinental® Spacer, 18mm Wide, Medium, 10°, 13mm	2	375.671	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 11mm	2
375.275	TransContinental® Spacer, 18mm Wide, Medium, 10°, 15mm	2	375.673	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 13mm	2
375.277	TransContinental® Spacer, 18mm Wide, Medium, 10°, 17mm	1	375.675	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 15mm	2
375.367	TransContinental® Spacer, 22mm Wide, Medium, 10, 7mm	2	375.677	TransContinental® Spacer, 18mm Wide, X-Large, 10°, 17mm	1
375.369	TransContinental® Spacer, 22mm Wide, Medium, 10°, 9mm	2	375.767	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 7mm	2
375.371	TransContinental® Spacer, 22mm Wide, Medium, 10°, 11mm	2	375.769	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 9mm	2
375.373	TransContinental® Spacer, 22mm Wide, Medium, 10°, 13mm	2	375.771	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 11mm	2
375.375	TransContinental® Spacer, 22mm Wide, Medium, 10°, 15mm	2	375.773	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 13mm	2

# TransContinental® 10° Lordotic Implant Set 975.924 (cont'd)

Part No.	Description	Qty	Part No.	Description
375.775	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 15mm	2	675.369	TransContinental® 22mm Trial, 10° Lordotic, 9mm
375.777	TransContinental® Spacer, 22mm Wide, X-Large, 10°, 17mm	1	675.371	TransContinental® 22mm Trial, 10° Lordotic, 11mm
375.867	TransContinental® Spacer, 18mm Wide, X-Small, 10°, 7mm	2	675.373	TransContinental® 22mm Trial, 10° Lordotic, 13mm
375.869	TransContinental® Spacer, 18mm Wide, X-Small, 10°, 9mm	2	675.375	TransContinental® 22mm Trial, 10° Lordotic, 15mm
375.871	TransContinental® Spacer, 18mm Wide, X-Small, 10°, 11mm	2	675.377	TransContinental® 22mm Trial, 10° Lordotic, 17mm
375.873	TransContinental® Spacer, 18mm Wide, X-Small, 10° 13mm	2	675.501	Insertion Sleeve
375.875	TransContinental® Spacer, 18mm	2	675.522 975.024	Insertion Sleeve, 22mm  TransContinental® 10° Lordotic
375.877	Wide, X-Small, 10°, 15mm  TransContinental® Spacer, 18mm  Wide, X-Small, 10°, 17mm	1	973.024	Implant Graphic Case
375.967	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 7mm	2		
375.969	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 9mm	2		
375.971	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 11mm	2		
375.973	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 13mm	2		
375.975	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 15mm	2		
375.977	TransContinental® Spacer, 22mm Wide, X-Small, 10°, 17mm	1		
675.065	TransContinental® 18mm Trial, 10° Lordotic, 5mm	1		
675.067	TransContinental® 18mm Trial, 10° Lordotic, 7mm	1		
675.069	TransContinental® 18mm Trial, 10° Lordotic, 9mm	1		
675.071	TransContinental® 18mm Trial, 10° Lordotic, 11mm	1		
675.073	TransContinental® 18mm Trial, 10° Lordotic, 13mm	1		
675.075	TransContinental® 18mm Trial, 10° Lordotic, 15mm	1		
675.077	TransContinental® 18mm Trial, 10° Lordotic, 17mm	1		
675.365	TransContinental® 22mm Trial, 10° Lordotic, 5mm	1		
675.367	TransContinental® 22mm Trial,	1		

10° Lordotic, 7mm

Qty

1

1

# TransContinental® Coronal Tapered Implant Set 975.925

Part No.	Description	Qty	Part No.	Description	Qty
375.092	TransContinental® Coronal Spacer, 0°, 18x30mm, 9mm	2	375.496	TransContinental® Coronal Spacer, 0°, 18x40mm, 17mm	1
375.093	TransContinental® Coronal Spacer, 0°, 18x30mm, 11mm	2	375.592	TransContinental® Coronal Spacer, 10°, 18x40mm, 9mm	2
375.094	TransContinental® Coronal Spacer, 0°, 18x30mm, 13mm	2	375.593	TransContinental® Coronal Spacer, 10°, 18x40mm, 11mm	2
375.095	TransContinental® Coronal Spacer, 0°, 18x30mm, 15mm	2	375.594	TransContinental® Coronal Spacer, 10°, 18x40mm, 13mm	2
375.096	TransContinental® Coronal Spacer, 0°, 18x30mm, 17mm	1	375.595	TransContinental® Coronal Spacer, 10°, 18x40mm, 15mm	2
375.192	TransContinental® Coronal Spacer, 10°, 18x30mm, 9mm	2	375.596	TransContinental® Coronal Spacer, 10°, 18x40mm, 17mm	1
375.193	TransContinental® Coronal Spacer, 10°, 18x30mm, 11mm	2	375.682	TransContinental® Coronal Spacer, 0°, 18x45mm, 9mm	2
375.194	TransContinental® Coronal Spacer, 10°, 18x30mm, 13mm	2	375.683	TransContinental® Coronal Spacer, 0°, 18x45mm, 11mm	2
375.195	TransContinental® Coronal Spacer, 10°, 18x30mm, 15mm	2	375.684	TransContinental® Coronal Spacer, 0°, 18x45mm, 13mm	2
375.196	TransContinental® Coronal Spacer, 10°, 18x30mm, 17mm	1	375.685	TransContinental® Coronal Spacer, 0°, 18x45mm, 15mm	2
375.292	TransContinental® Coronal Spacer, 0°, 18x35mm, 9mm	2	375.686	TransContinental® Coronal Spacer, 0°, 18x45mm, 17mm	1
375.293	TransContinental® Coronal Spacer, 0°, 18x35mm, 11mm	2	375.692	TransContinental® Coronal Spacer, 0°, 18x50mm, 9mm	2
375.294	TransContinental® Coronal Spacer, 0°, 18x35mm, 13mm	2	375.693	TransContinental® Coronal Spacer, 0°, 18x50mm, 11mm	2
375.295	TransContinental® Coronal Spacer, 0°, 18x35mm, 15mm	2	375.694	TransContinental® Coronal Spacer, 0°, 18x50mm, 13mm	2
375.296	TransContinental® Coronal Spacer, 0°, 18x35mm, 17mm	1	375.695	TransContinental® Coronal Spacer, 0°, 18x50mm, 15mm	2
375.392	TransContinental® Coronal Spacer, 10°, 18x35mm, 9mm	2	375.696	TransContinental® Coronal Spacer, 0°, 18x50mm, 17mm	1
375.393	TransContinental® Coronal Spacer, 10°, 18x35mm, 11mm	2	375.782	TransContinental® Coronal Spacer, 10°, 18x45mm, 9mm	2
375.394	TransContinental® Coronal Spacer, 10°, 18x35mm, 13mm	2	375.783	TransContinental® Coronal Spacer, 10°, 18x45mm, 11mm	2
375.395	TransContinental® Coronal Spacer, 10°, 18x35mm, 15mm	2	375.784	TransContinental® Coronal Spacer, 10°, 18x45mm, 13mm	2
375.396	TransContinental® Coronal Spacer, 10°, 18x35mm, 17mm	1	375.785	TransContinental® Coronal Spacer, 10°, 18x45mm, 15mm	2
375.492	TransContinental® Coronal Spacer, 0°, 18x40mm, 9mm	2	375.786	TransContinental® Coronal Spacer, 10°, 18x45mm, 17mm	1
375.493	TransContinental® Coronal Spacer, 0°, 18x40mm, 11mm	2	375.792	TransContinental® Coronal Spacer, 10°, 18x50mm, 9mm	2
375.494	TransContinental® Coronal Spacer, 0°, 18x40mm, 13mm	2	375.793	TransContinental® Coronal Spacer, 10°, 18x50mm, 11mm	2
375.495	TransContinental® Coronal Spacer, 0°, 18x40mm, 15mm	2	375.794	TransContinental® Coronal Spacer, 10°, 18x50mm, 13mm	2

# TransContinental® Coronal Tapered Implant Set 975.925 (cont'd)

Part No.	Description	Qty	Part No.	Description	Qty
375.795	TransContinental® Coronal Spacer, 10°, 18x50mm, 15mm	2	675.031	Trial Leading – 4° Coronal x 0° Lordotic, 11mm	1
375.796	TransContinental® Coronal Spacer, 10°, 18x50mm, 17mm	1	675.033	Trial Leading – 4° Coronal x 0° Lordotic, 13mm	1
375.882	TransContinental® Coronal Spacer, 0°, 18x55mm, 9mm	2	675.035	Trial Leading – 4° Coronal x 0° Lordotic, 15mm	1
375.883	TransContinental® Coronal Spacer, 0°, 18x55mm, 11mm	2	675.037	Trial Leading – 4° Coronal x 0° Lordotic, 17mm	1
375.884	TransContinental® Coronal Spacer, 0°, 18x55mm, 13mm	2	675.227	Trial Trailing – 4° Coronal x 0° Lodotic, 7mm	1
375.885	TransContinental® Coronal Spacer, 0°, 18x55mm, 15mm	2	675.229	Trial Trailing – 4° Coronal x 0° Lodotic, 9mm	1
375.886	TransContinental® Coronal Spacer, 0°, 18x55mm, 17mm	1	675.231	Trial Trailing – 4° Coronal x 0° Lodotic, 11mm	1
375.892	TransContinental® Coronal Spacer, 0°, 18x60mm, 9mm	2	675.233	Trial Trailing – 4° Coronal x 0° Lodotic, 13mm	1
375.893	TransContinental® Coronal Spacer, 0°, 18x60mm, 11mm	2	675.235	Trial Trailing - 4° Coronal x 0° Lodotic, 15mm	1
375.894	TransContinental® Coronal Spacer, 0°, 18x60mm, 13mm	2	675.237	Trial Trailing - 4° Coronal x 0° Lodotic, 17mm	1
375.895	TransContinental® Coronal Spacer, 0°, 18x60mm, 15mm	2	675.427	Trial Leading – 4° Coronal x 10° Lordotic, 7mm	1
375.896	TransContinental® Coronal Spacer, 0°, 18x60mm, 17mm	1	675.429	Trial Leading – 4° Coronal x 10° Lordotic, 9mm	1
375.982	TransContinental® Coronal Spacer, 10°, 18x55mm, 9mm	2	675.431	Trial Leading – 4° Coronal x 10° Lordotic, 11mm	1
375.983	TransContinental® Coronal Spacer, 10°, 18x55mm, 11mm	2	675.433	Trial Leading – 4° Coronal x 10° Lordotic, 13mm	1
375.984	TransContinental® Coronal Spacer, 10°, 18x55mm, 13mm	2	675.435	Trial Leading – 4° Coronal x 10° Lordotic, 15mm	1
375.985	TransContinental® Coronal Spacer, 10°, 18x55mm, 15mm	2	675.437	Trial Leading – 4° Coronal x 10° Lordotic, 17mm	1
375.986	TransContinental® Coronal Spacer, 10°, 18x55mm, 17mm	1	675.827	Trial Trailing – 4° Coronal x 10° Lordotic, 7mm	1
375.992	TransContinental® Coronal Spacer, 10°, 18x60mm, 9mm	2	675.829	Trial Trailing – 4° Coronal x 10° Lordotic, 9mm	1
375.993	TransContinental® Coronal Spacer, 10°, 18x60mm, 11mm	2	675.831	Trial Trailing – 4° Coronal x 10° Lordotic, 11mm	1
375.994	TransContinental® Coronal Spacer, 10°, 18x60mm, 13mm	2	675.833	Trial Trailing – 4° Coronal x 10° Lordotic, 13mm	1
375.995	TransContinental® Coronal Spacer, 10°, 18x60mm, 15mm	2	675.835	Trial Trailing – 4° Coronal x 10° Lordotic, 15mm	1
375.996	TransContinental® Coronal Spacer, 10°, 18x60mm, 17mm	1	675.837	Trial Trailing – 4° Coronal x 10° Lordotic, 17mm	1
675.027	Trial Leading – 4° Coronal x 0° Lordotic, 7mm	1	675.940	Holder Inserter	2
675.029	Trial Leading – 4° Coronal x 0°	1	675.960	Coronal Insertion Sleeve, 18mm	2
	Lordotic, 9mm		975.025	TransContinental® Coronal Tapered Implant Graphic Case	

# IMPORTANT INFORMATION ON MARS™ (Minimal Access Retractor System)

# **DESCRIPTION**

MARS™ (Minimal Access Retractor System) is a comprehensive retractor, ports and instrument system that provides efficient access to posterior lumbar spine. MARS™ and MARS™3V consists of a retractor frame, blades, disposable ports, silicone sleeves, light cables and associated manual surgical instruments. The blades and ports are available in several designs to accommodate individual patient anatomy.

The MARS™ and MARS™3V instruments are made from aluminum and stainless steel as specified in ASTM B221-02 and F899-02. The ports are made from radiolucent polymer (PEEK) as specified in ASTM F2026.

# **CLEANING**

Cleaning instructions by hand, when properly carried out, causes less damage than mechanical cleaning. When cleaning instruments by hand, the following should be observed:

- 1. Clear any corners or recesses of all debris. (Note: extra care should be taken to clean out any cannulated areas by using an appropriate cleaning stylet and rinsing immediately.)
- 2. Remove all traces of blood and other such residues immediately. Do not allow these to dry.
- 3. The instruments should be submerged (if applicable) and cleaned with a commercially available manual cleaner (i.e. Instraclean from Calgon or Medline High Suds Detergent) prepared according to the manufacturer's recommendation.
- 4. A soft nylon bristled brush is then used to manually clean the devices while immersed in the cleaning solution. Never use steel brushes or abrasive pads, as these rupture the passive layer of the instrument surface which can lead to corrosion.
- 5. The instruments should be thoroughly rinsed after cleaning. Distilled water should be used.
- 6. Dry instruments immediately after cleaning.

These devices are supplied NONSTERILE. Sterilization is recommended as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Gravity Displacement (Wrapped)	132°C (270°F)	25 Minutes	45 Minutes
Steam	Pre-vacuum (Wrapped)	132°C (270°F)	15 Minutes	30 Minutes

Cycles should be performed on tray with devices opened for maximum

These parameters are validated to sterilize only these instruments. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

The following information is provided by LumitexMD, Inc. for MARS™ Light cable distributed by Globus Medical, Inc.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

# INFORMATION FOR USE FOR THE MARS™ LIGHT CABLE

# DESCRIPTION

The MARS™ Light Cable is a sterile, single use, latex free, plastic fiber optic device intended to bring cool area lighting into deep surgical sites. The MARS™ Light Cable is intended for use with a 300 watt xenon illuminator, using a 4mm fiber optic cable with a female ACMI connector. For best results use a Globus cable by LumitexMD.

# INDICATIONS FOR USE

The MARS™ Light Cable is intended for the illumination of surgical procedures, particularly where deep cavities or adjacent tissues limit outside light in the surgical field. It is designed for use in less invasive spinal surgery.

# CONTRAINDICATIONS

The MARS™ Light Cable presents no contraindication. However, the user should be familiar with the use of light sources and cables and should take precautions accordingly.

# **WARNINGS**

The MARS™ Light Cable is designed for use with 300 watt xenon illuminators, using a 4mm fiber optic cable. Do not use light sources rated higher than 300 watts, or cables with fiber optic bundles of more than 4mm diameter. Use of higher watt sources or larger diameter cables could result in overheating; causing product failure and patient injury.

Should the MARS™ Light Cable become cut, collect fluid inside, appear broken or damaged in any manner, it should be replaced to minimize risk to the patient.

Do not operate the light source and cable without the MARS™ Light Cable attached. Without the MARS™ Light Cable, the output from the fiberoptic cable is extremely bright, hot and may cause burns, ignite drapes/gowns, or temporarily blind vision.

# **PRECAUTIONS**

Light sources vary widely in emission of visible and infrared energy. As a precautionary measure, we recommend occasionally monitoring connector temperature during first time use with a new light source or lamp; thereafter if needed. As is common with fiber optic equipment, metal portion of connector can become hot to the touch. Use plastic grip as handle. Do not place the metal ring portion of connector directly on the patient's skin.

Because light energy can be absorbed as heat, the entire lit portion (distal end) of the MARS™ Light Cable should not be continuously embedded (i.e. lit surface should not be completely buried) in tissue and held fixed for more than a few minutes at one time.

Each MARS™ Light Cable package contains one MARS™ Light Cable assembly with an integrated adhesive strip and two double-sided adhesive strips. Each adhesive strip includes two paper release liners. Prior to closing the surgical site, all components must be accounted for.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

# IMPORTANT INFORMATION ON MARS™ (Minimal Access Retractor System)

# **DIRECTIONS FOR USE**

Attach the MARS™ Light Cable to the Globus Medical MARS™ retractor using the integrated stainless steel clip located on the back of each MARS™ Light Cable.

The MARS™ Light Cable connects to a light source used for head lamps or endoscopes. A fiber optic cable attaches the light source and MARS™ Light Cable. Make sure the MARS<sup>™</sup> Light Cable connector is securely attached to the cable. The cable should be in good repair with clean optics. Dirty optics or cables in need of repair can cause excessive heat at the connectors.

Turning down overhead lighting may improve visualization within the surgical site.

Body fluids or debris collecting on the surface of the MARS™ Light Cable may be irrigated or wiped away.

Sterile unless package is opened or damaged. Do not use if package is opened or damaged.

# LIMITED WARRANTY

LumitexMD warrants the material conformity of the MARS™ Light Cable to specifications in the product labeling until the earlier of 12 months from shipment to customer or the expiration date of the product, and will repair or replace at LumitexMD option and expense any LumitexMD product that does not meet specifications in all material respects. LUMITEXMD LIABILITY TO CUSTOMER, USER, OR PATIENT IS EXPRESSLY LIMITED TO REPAIR OR REPLACEMENT. LumitexMD expressly disclaims all other warranties, express or implied, including, without limitation, merchantability or fitness for a particular purpose. Please direct any inquires to Globus Medical.

Distributed by: Globus Medical, Inc. Valley Forge Business Center 2560 8443 Dow Circle General Armisted Avenue Audubon, PA 19403 Tel: 610-415-9000 Fax: 610-415-9144


Manufactured by: LumitexMD, Inc. Tel: 800-969-5483 www.lumitexmd.com medical@lumitex.com

Authorized EC Representative: Medical Device Safety Service GmbH Schiffgraben 41 Strongsville, OH 44136 USA D-30175 Hannover, Germany









# IMPORTANT INFORMATION ON THE PATRIOT® LUMBAR SPACER SYSTEM

# **DESCRIPTION**

The PATRIOT® Spacers (CONSTITUTION® PLIF, SIGNATURE® TLIF, CONTINENTAL® ALIF, TransContinental® and TransContinental® M Spacers) are lumbar interbody fusion devices used to provide structural stability in skeletally mature individuals following discectomy. Each of the PATRIOT® spacers provides a different shape to accommodate various surgical approaches to the lumbar spine. The CONSTITUTION® PLIF Spacer is inserted using a posterior approach. The SIGNATURE® TLIF Spacer is inserted using a transforaminal approach. The CONTINENTAL® ALIF Spacer is inserted using an anterior approach. The TransContinental® and TransContinental® M Spacer are inserted using an anterior or lateral approach. The devices are available in various heights and geometric options to fit the anatomical needs of a wide variety of patients. These spacers are to be filled with autogenous bone graft material. Protrusions on the superior and inferior surfaces of each device grip the endplates of the adjacent vertebrae to resist expulsion.

PATRIOT® Spacers are made from PEEK radiolucent polymer, with titanium alloy or tantalum markers, as specified in ASTM F2026, F136, F1295, and F560. The TransContinental® M Spacer also includes an integrated titanium alloy nut, as specified in ASTM F136, F1295.

# **INDICATIONS**

PATRIOT® Spacers (CONSTITUTION® PLIF, SIGNATURE® TLIF, CONTINENTAL® ALIF, TransContinental® and TransContinental® M Spacers) are interbody fusion devices intended for use in patients with degenerative disc disease (DDD) at one or two contiguous levels of the lumbosacral spine (L2-S1). DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had at least six (6) months of non-operative treatment. In addition, these patients may have up to Grade 1 spondylolisthesis or retrolisthesis at the involved level(s).

PATRIOT® Spacers are to be filled with autogenous bone graft material. These devices are intended to be used with supplemental fixation.

# **WARNINGS**

One of the potential risks identified with this system is death. Other potential risks which may require additional surgery, include:

- device component fracture,
- loss of fixation,
- non-union,
- fracture of the vertebrae,
- neurological injury, and
- vascular or visceral injury.

Interbody fusion devices for the treatment of degenerative conditions are designed to withstand both full load bearing and the loads associated with long-term use which could result from the presence of non-union or delayed union.

Certain degenerative diseases or underlying physiological conditions such as diabetes, rheumatoid arthritis, or osteoporosis may alter the healing process, thereby increasing the risk of implant breakage or spinal fracture.

Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery.

Components of this system should not be used with components of any other system or manufacturer.

The components of this system are manufactured from PEEK radiolucent polymer, titanium alloy, and tantalum. Mixing of stainless steel implant components with different materials is not recommended for metallurgical, mechanical and functional reasons.

# **PRECAUTIONS**

The implantation of intervertebral fusion devices should be performed only by experienced spinal surgeons with specific training in the use of this system because this is a technically demanding procedure presenting a risk of serious injury to the patient. Preoperative planning and patient anatomy should be considered when selecting implant size.

Surgical implants must never be reused. An explanted implant must never be reimplanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which could lead to breakage.

Adequately instruct the patient. Mental or physical impairment which compromises or prevents a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

The PATRIOT® Spacers have not been evaluated for safety and compatibility in the MR environment. The PATRIOT® Spacers have not been tested for heating or migration in the MR environment.

For optimal implant performance, when using the PATRIOT® Spacers, the physicians/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of this system.

# CONTRAINDICATIONS

Use of PATRIOT® Spacer(s) is contraindicated in patients with the following

- 1. Active systemic infection, infection localized to the site of the proposed implantation, or when the patient has demonstrated allergy or foreign body sensitivity to any of the implant materials.
- 2. Prior fusion at the level(s) to be treated.
- 3. Severe osteoporosis, which may prevent adequate fixation
- 4. Conditions that may place excessive stresses on bone and implants, such as severe obesity or degenerative diseases, are relative contraindications. The decision whether to use these devices in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
- 5. Patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stresses on the implant during bony healing and may be at a higher risk of implant failure.
- 6. Any condition not described in the indications for use.

# **CONTACT INFORMATION**

Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

# **STERILIZATION**

The PATRIOT® Spacer implants and instruments have been validated to assure a Sterility Assurance Level (SAL) of 10-6. The use of an FDA cleared wrap is recommended, per the Association for the Advancement of Medical Instrumentation (AAMI) ST79, Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

# Implants:

These devices are supplied NONSTERILE. Sterilization is recommended as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Gravity Displacement (Wrapped)	132°C (270°F)	10 minutes	15 minutes
Steam	Pre-vacuum (Wrapped)	132°C (270°F)	4 minutes	15 minutes

# Instruments:

These instruments are supplied NONSTERILE. Sterilization is recommended as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Gravity Displacement (Wrapped)	132°C (270°F)	25 minutes	15 minutes
Steam	Pre-vacuum (Wrapped)	132°C (270°F)	15 minutes	20 minutes

These parameters are validated to sterilize only this device. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.


CAUTION: Federal (USA) Law Restricts this Device to Sale by or on the order of a Physician.



Notes				

# Notes

Notes				





Globus Medical Valley Forge Business Center 2560 General Armistead Avenue Audubon, PA 19403 www.globusmedical.com

Customer Service:

Phone 1-866-GLOBUS1 (or 1-866-456-2871) 1-866-GLOBUS3 (or 1-866-456-2873) Fax

© 2012 Globus Medical. All rights reserved. Patents pending. Life moves us is a registered trademark of Globus Medical. MARS is a trademark of Globus Medical. CONTINENTAL, PATRIOT, REVERE, REVOLVE and TransContinental are registered trademarks of Globus Medical. Please refer to package insert for description, indications, contraindications, warnings, precautions and other important information.