


Speedband Superview Super 7™

Boston
Scientific

Multiple Band Ligator




Ergonomic Handle

Single String Deployment

Speedband Superview Super 7™

Multiple Band Ligator

Banding has been shown to be effective in the treatment of esophageal variceal bleeding and anorectal hemorrhoids. The Speedband Superview Super 7 Multiple Band Ligator is designed to provide clear visibility and enhanced band performance in a user-friendly system.



CLEAR LIGATING UNIT WITH REAR-MOUNTED BANDS

- Provides excellent visibility for rapid visual navigation of anatomical landmarks and pathology.

ALERT BAND

- Natural colored band is designed to alert the operator when only one band remains.



SINGLE STRING DEPLOYMENT

- Engineered to provide excellent suction for tissue capture, band placement and band retention.

SQUARE BAND DESIGN

- Molded bands are designed for consistent shape and reliable performance.
- Latex-free material is designed to reduce the risk of allergic reaction in sensitive patients.

REDESIGNED ERGONOMIC HANDLE

- User-friendly set-up.

Engineered for Operator Comfort

ORDERING INFORMATION

LATEX

Order Number	Product Description	Number of Bands	Units
M00542251	Multiple Band Ligator	7	Box 2
M00542253	Multiple Band Ligator	7	Box 4

Boston Scientific

Advancing science for life™

Boston Scientific Corporation
100 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com/endo-resources

Ordering Information
1.888.272.1001

©2017 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-153507-AB August 2017

Registered trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.