

Segura[™] Baskets

Stone Retrieval Baskets

Boston
Scientific

Permaform[™] Wire
Technology gives
Segura Baskets an
added advantage



Segura™ Baskets

Permaform Flatwire Construction

- Designed to maintain basket symmetry during stone capture even when difficult manipulations are required.
- Produces exceptional wire strength for secure grasp and reliable retrieval.
- Permits manufacture of versatile three-wire basket design.

Exclusive 13mm Three-Wire Basket Design

- Large 120° opening combined with 13mm basket O.D. is ideal for retrieval of large stones such as those encountered in the renal calyces.
- Opening force is equivalent to that of full-sized Segura Baskets.

Positive Action Detachable Handle

- Ergonomic design provides excellent feedback for precise single hand control.

Ordering Information

Segura Hemisphere™ Stone Retrieval Baskets (rounded edge wire)

Order Number	Description	Sheath O.D. (Fr)	Working Length (cm)	Basket O.D. (mm)	Sheath Material
M0063801000	4 Wire Basket, "O" Tip	3.0	120	16	Polyimide
M0063801020	4 Wire Basket, "O" Tip	3.0	90	16	Polyimide
M0063801060	4 Wire Basket, "O" Tip	2.4	90	16	Polyimide
M0063801070	4 Wire Basket, "O" Tip	2.4	120	16	Polyimide
M0063801100	4 Wire Basket, "O" Tip	3.0	120	20	Polyimide
M0063801120	4 Wire Basket, "O" Tip	3.0	90	20	Polyimide
M0063801160	4 Wire Basket, "O" Tip	4.5	90	16	PTFE
M0063801180	4 Wire Basket, "O" Tip	4.5	90	20	PTFE
M0063801190	4 Wire Basket, 5cm Filiform Tip	4.5	90	20	PTFE

Segura Baskets (flatwire)

Order Number	Description	Sheath O.D. (Fr)	Working Length (cm)	Basket O.D. (mm)
M0063003110	3 Wire Basket, "O" Tip	3.0	120	13
M0063001000	6 Wire Basket, "O" Tip	4.5	90	20
M0063001010	6 Wire Basket, 5cm Filiform Tip	4.5	90	20

Segura-Dretler™ Laser Baskets* (flatwire)

Order Number	Description	Recommended Fiber Core (µ)	Sheath O.D. (Fr)	Working Length (cm)	Basket O.D. (mm)
M0063201030	Segura-Dretler, 4 Wire Basket, "O" Tip	320	4.5	90	20
M0063201040	Segura-Dretler, 4 Wire Basket, "O" Tip	320	3.0	90	16

*Not for use with Holmium laser fibers.

Over Guidewire Segura Basket (flatwire)

Order Number	Description	Sheath O.D. (Fr)	Working Length (cm)	Basket O.D. (mm)
M0063201020	Over Guidewire, 4 Wire Basket, "O" Tip	7.0	70	20

Supplied With: .038" Stiff Bodied Stainless Steel Guidewire

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician. Refer to the Directions for Use provided with this product for complete instructions, warnings, and precautions prior to using this product.

PERMAFORM WIRE TECHNOLOGY

Permaform Wire Segura Baskets, made possible by proprietary design and manufacturing techniques, establish a standard of performance for flatwire baskets.

Shape Stability

Permaform Wire construction is designed to minimize the tendency of traditional flatwire baskets to invert during difficult stone manipulations.

Wire Strength

Permaform Wire construction is designed to provide excellent strength for reliable opening force and secure grasp.

Segura Hemisphere Basket Configurations

Unique Rounded Edge Wire is designed to reduce ureteral trauma during stone manipulation.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.800.225.3226

Copyright © 2003 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

DVU860 5M 4/03-4/05