

Ureteroscopy Accessories

Leak-resistant accessories designed to make your cases flow smoothly

SAPS™ Single Action Pumping System

A pumping system designed for physicians or assistants to effectively control irrigation and clear the field of view with a simple one-handed press of the syringe.

- One-way valves with automatic refill of the syringe for convenient, immediate access to irrigation fluid



SAPS™ CF Continuous Flow version pictured



UroLok™ II Scope Adaptor

Our UroLok Scope Adaptor prevents backflow and secures slippery, coated devices. It connects directly to drainage bags and is designed to keep these connections aside and out of the way.

Gateway™ Advantage Y-Adaptor

A scope adaptor designed with a leak-resistant seal and two ports to smoothly and securely introduce multiple instruments into the scope.

- Main port allows passage of devices up to 8F
- Blue luer seal accepts devices up to 3F



Ordering Information

SAPS™ CF Continuous Flow Single Action Pumping System

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|---------------------------------|------|
| M0067201011 | 08714729045076 | Single Action Pumping System CF | bx 5 |

System Includes: 10mL Vacuum Syringe, Continuous Flow Adjustable Length Tubing System with Dual-Check Valve and Roller Clamp

UroLok™ II Scope Adaptor

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|-------------------------|-------|
| M0067301401 | 08714729066569 | UroLok II Scope Adaptor | bx 10 |

Encore™ 26 Inflator

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|---|------|
| M0067101131 | 08714729761815 | 20mL Encore 26 Inflator with Pressure Gauge | bx 5 |
| M0067101140 | 08714729755814 | 20mL Encore 26 Inflator with Pressure Gauge | ea |

Ureteral Catheter Connectors

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|---|------|
| M0067401100 | 08714729010982 | Connection for 3 through 6F Catheters with Luer Lock Fitting | ea |
| M0067401110 | 08714729010975 | Connection for 7 through 10F Catheters with Luer Lock Fitting | ea |

SAPS™ Single Action Pumping System

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|------------------------------|------|
| M0067201001 | 08714729045069 | Single Action Pumping System | bx 5 |

System Includes: 10mL Vacuum Syringe, Check Valve with Adjustable Length Connecting Tubes

Gateway™ Advantage Y-Adapter

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|-----------------------------|-------|
| M0067301501 | 08714729845331 | Gateway Advantage Y-Adapter | bx 10 |

Connecting Tubes for Drainage Bag

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|---|------|
| M0067401000 | 08714729010944 | 30cm Connecting Tube for 3 through 6F Catheters | ea |
| M0067401010 | 08714729012832 | 30cm Connecting Tube for 7 through 10F Catheters | ea |
| M0067401020 | 08714729010951 | 30cm "Y" Connecting Tube for two 3 through 6F Catheters | ea |

Drainage Bag Connectors

| Order Number | GTIN Number | Description | Unit |
|--------------|----------------|--|------|
| M0067401200 | 08714729010999 | Connection for Female Luer Lock Fitting to Drainage Bag | ea |
| M0067401220 | 08714729011002 | Connection for 3 through 6F Ureteral Catheters to Drainage Bag | ea |

Applicable Medicare C-codes: There are no applicable Medicare C-codes for the above products.

To order, contact your sales representative.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

©2017 Boston Scientific Corporation or its affiliates. All rights reserved.

URO-446715-AA JAN 2017